SIMPLEE®

All you need. No more. But no less.

HELIX TEST

GB Instructions for use

Simplee Helix Test is a sterilization indicator and batch control system designed for testing Class B and S steam autoclaves. PCD complies to NEN-EN 867-5. Indicators complies to ISO 11140-1 class 2.

- Ensure the Helix Tube is dry and not warm. Unscrew the cap from the Helix Tube.
- Take an indicator strip out of the package (image 1) and fold it in the middle with the text on the inside (image 2).
- Insert the middle fold part of the indicator strip first into the closure cap as far as it will go into the cap (image 3).
- Screw the closure cap back into the Helix Tube (image 4) and place in the autoclave.
- Run a sterilization process as indicated on the strip: 134°C - 3.5 min or 121°C - 15 min.
- At the end of the sterilization process, remove the tube from the autoclave. Be careful, the tube may be hot!
- Unscrew the cap and check the indicator strip. If the indicator shows green (image 5) it can be assumed that

the autoclave is functioning properly and sterilizes successfully.

- The indicator strip is provided with an adhesive layer so that the strips with additional information about the process and the load can be easily archive.

Note

Test Kit should be renewed after 250 batches.

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Lifco Dental AB; Verkämästaregatan 1 · 745 39 Enköping/Sweden