SIMPLEE®

All you need. No more. But no less.

RECIPFILE V

(GB)

Specifications

Model

25#, 40#, 50#, 25-50#

Length

21mm, 25mm, 31mm

Material

NiTi alloy

Advantages

- Flexibility
- excellent cutting efficiency
- Reduce the risk of contamination
- Simplicity, one file system

Recommend Speed and Torque

Model	Speed (RPM)	Taper
25#	250 - 350	.08
40#	250 - 350	.06
50#	250 - 350	.05

Motion

Reciprocating movement

Package

Blister package 6 pcs/blister

Sequence

- To achieve a straight line access to the root canal entrance.
- Select the correct Simplee V files instruments
- Check the Endo motor has been selected
- Introduce the Simplee V files into the root canal.
- Move the instrument in a slow in-and-out pecking motion. The
- amplitude of the in-and-out movements should not exceed 3 mm.
- Only light pressure should be applied. 1 in-and-out movement=1 peck.
- The instrument will advance easily in the canal.
- After 3 pecks, remove the instrument from the canal. Clean the debris.
- Irrigate the canal.
- Make sure the canal is free to approx 3mm beyond the prepared
- canal section with an ISO size 10 file.
- In this way, continue with the Simplee V files instrument until approx 2/3
- of the working length has been reached.
- When using an 25#
- determine the working length by using an ISO size 10# files. When
- using an 40#, 50#, the working length should be re-checked with an
- apex locator.
- Continue with the Simplee V files instrument until working length has
- been reached.
- As soon as full working length has been reached, withdraw the
- instrument from the root canal.

Selecting the correct Simplee V files

- If the canal is partially or completely invisible on the radiograph:
- the canal is considered narrow, use size 25#
- If the canal is completely visible on the radiograph: take an ISO
- 30# hand file, insert it passively to working length. If it reached
- working length, the canal is considered large, use 50#.
- If an ISO 30 hand files does not goes passively to working length,
- try an ISO 20# hand file. If this goes passively to working length, the
- canal is considered medium, use 40#.

Technique tips

- Simplee V files instruments can be used in a lateral brushing motion to enable preparation of irregular shaped canals to enlarge the root canal entrance. The lateral brushing motion can also help to achieve an easier advancement of the instrument.
- Clean Simplee V files instrument flute after 3 pecks.
- Use an ISO size 10# file to check the situation of the canal, confirm the canal is not blocked after 3 pecks with Simplee V files.
- Irrigate the root canal according to the appropriate irrigation protocol.
- Never apply pressure if resistance is encountered, instead clean the files and check the canal with an ISO file size 10#.
- Remove the instrument from the canal as soon as working length has been reached. Working with engine files too long in one spot will cause canal transportation.
- Pull the Simplee V files out of the canal after 3 pecks or when resistance is encountered.

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