# SIMPLEE®

All you need. No more. But no less.

## RECIPFILE VB

(GB)

### Specifications

Model: 25#, 40#, 50#, 25-50#

### Material

NiTi alloy with heat treatment to increase the resistance to cyclic fatigue.

### Advantages

- Super flexibility
- excellent cutting efficiency
- Reduce the risk of contamination
- one file system, time-saving and easy to use
- Minimized the fracture of instrument
- Re-treatment

### Specifications

Model	Length(mm)	taper
25#	21,25,31	.08
40#	21,25,31	.06
50#	21,25,31	.05

### Motion

Reciprocating movement

### Package

Blister package 6pcs/blister

### Sequence

- To achieve a straight line access to the root canal entrance.
- Irrigating the root canal.
- Select the correct Simplee RecipFile VB instruments
- Check the Endo motor has been selected
- Introduce the Simplee RecipFile VB into the root canal.
- Move the instrument in a slow in-and-out pecking motion. The amplitude of the in-and-out movements should not exceed 3mm.
- Only light pressure should be applied. 1 in-and-out movement = 1 peck. The instrument will advance easily in the canal.
- After 3 pecks, remove the instrument from the canal. Clean the debris.
- Irrigate the canal.
- Make sure the canal is free to approx 3mm beyond the prepared canal section with an ISO size 10 file.
- Continue with the Simplee RecipFile VB instrument until working length has been reached.
- As soon as full working length has been reached, withdraw the instrument from the root canal.

### Selecting the correct Simplee RecipFile VB

- If the canal is partially or completely invisible on the radiograph: the canal is considered narrow, use size 25#
- If the canal is completely visible on the radiograph: take an ISO 30# hand file, insert it passively to working length. If it reached working length, the canal is considered large, use 50#.
- If an ISO 30 hand files does not goes passively to working length, try an ISO 20# hand file. If this goes passively to working length, the canal is considered medium, use 40#.

### Technique tips

- Pull the instrument out of the canal after 3 pecks or when resistance is encountered.
- Never apply pressure if the instrument does not advance into the canal.
- Don't forget to clean and inspect the flutes during preparation.
- Remove the file from the canal as soon as working length has been reached.
- Simplee RecipFile VB can be used in lateral brushing motion to enable preparation of irregularly shaped canals.
- It is not recommended to use Simplee RecipFile VB nor continuous rotary system if the root canal has an abrupt curvature in the apical section.

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