SIMPLEE®

All you need. No more. But no less.

W+ SYSTEM

(GB)

Directions for use

FOR DENTAL USE ONLY
RECIPROCATING MOTION

RAW MATERIAL

The whole cutting area of W+ system files are made of
Nickel Titanium alloy.

SPECIFICATIONS

- Simplee W+ Tunnel-Feile 15/02
- Simplee W+ Shaping-Feilen

Small W+	20/07
Primary W+	25/07
Medium W+	35/06
Large W+	45/05

FUNCTION

- Simplee W+ Tunnel files are used in the purpose of building a smooth tunnel during preparation by mechanically powered W+ files system.
- Simplee W+ Shaping files are used aiming at shaping and cleaning the root canal by mechanically powered W+ files system

WARNINGS

- Please confirm that patient has no allergy to Nickel Titanium alloy and other medical reagent before the treatment.
- These instruments are intended for use only in dental root canal treatment, users must follow a correct and reasonable way under medical environments.
- Please use a rubber dam system to ensure patients' safety when deliver endodontic treatment.
- Simplee W+ system are all disposable instruments, to avoid cross contamination and ensure it can reach a optimum performance, please abandon it after a single use.

PRECAUTIONS

- Safety equipment (gloves, glasses, rubber dam) is required for patients' and self-protection.
- These instruments can only be used after sterilization, please using these files under medical environment.
- Simplee W+ system are all heat-treated instruments and showing a greater flexibility performance than Simplee W system, it's a normal phenomenon that it might appear a little bent but it won't bring any disturbance for use.
- Reciprocating motion is recommended to match the Simplee W+ system.
- Watch out for the disposable warnings, in order to get a ideal performance of these instruments please not to try to use more than once.
- Please use Simplee W+ Tunnel file and Simplee W+ Shaping file in a appropriate case and in a correct way.
- Before using Simplee W+ Tunnel file, scouting by hand files is required and necessary.
- Remove debris out frequently from screwing-in problems
- Normally shaping the canal with EDTA for entering the canal smoothly and irrigate sufficiently after removing debris created by shaping the root canal.
- If you come across a situation that it's difficult to enter the deeper part of the canal, please stop and check both morphology of the canal and whether there are signs of deformation of files or not, then confirm the orifice is expanded wide enough for subsequent use of shaping files.

OPERATING STEPS

- Preoperative radiographic evaluation**
Observe the configuration of canal, then confirm the quantity, width, length, curvature, and choose a most suitable file system to deliver a treatment.
- Access**
Prepare a straight access to canal orifice.
- Scouting**
Explore the canal up to size #10 hand use file with EDTA.
- Smooth path done by Simplee W+ Tunnel file**
Create a smooth path and confirming the whole working length can be totally reached after the orifice is well expanded in order to be well prepared for subsequent use of shaping files, as a excellent connection between hand use files and rotary files.
- Shaping with Simplee W+ Shaping file**
On the basis of created glide path, shape the canal with Simplee W+ shaping file in the condition of EDTA. Normally choose Primary file (25/07) as the initial file to start, changing the most optimal file subject to morphology of the canal.

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