SIMPLEE®

All you need. No more. But no less.

INSTRUCTIONS OF USE & STERILIZATION

Warnings

Instruments are not sterile upon receipt and must be sterilized before use in accordance with the following instruction. Instruments can be damaged by alkaline and acidic detergents. Always use a pH (7) neutral detergent or solution. Temperature do not exceed 134°C and remove all packaging prior to cleaning and sterilization, contaminated instruments should be handled with gloves and eye protection should be used.

Limitations

Repeated processing has minimal effect on these instruments. End of life is normally determined by wear and damage due to use. Never clean color-etched probes or black composite instruments in ultrasonic baths.

Preparation for Cleaning

No specific requirement, Disassemble not required. Rinsing: Immediately after surgery rinse instruments under warm (not hot) water. It may be helpful to use a nylon toothbrush to rinse the lock boxes and joints of the instruments. Be sure to remove all blood, body fluids, and tissue. Cleaning: If you do not clean your instrument immediately after rinsing, instruments should be submerged in a solution of water and neutral pH (7) detergent. They should never be placed in saline solution, as it may cause corrosion and eventually irreversible damage to the instrument.

Ultrasonic Cleaning

For micro and delicate instruments, use manual cleaning. Instruments should be processed in a cleaner for the full recommended cycle time – usually 5 to 10 minutes. Place instruments in open position into the ultrasonic cleaner. Make sure that "Sharp" instruments do not touch other instruments. All instruments have to be fully submerged. Do not place dissimilar metals (stainless, copper, chrome, plated, etc.) in the same cleaning cycle. Change solution freq.

Automated reprocessing

Simplee instruments are suitable for machine reprocessing in the thermal disinfectant (WD). Please follow manufacturer's instructions for the respective thermal disinfectant (WD).

After Cleaning

If instruments are to be stored, let them, air dry and store them in a clean and dry environment.

Thermal disinfection

Instead of manual cleaning and disinfection, the instruments can be reprocessed in an automated washer-disinfectant (WD). Please follow the specific manufacturer's instructions for the washer-disinfectant.

Autoclaving

If instruments are to be reused or autoclaved maintain lubricate of all instruments which have any "metal to metal" action such as scissors, Hemostats, needle holder, self-retaining retractors, etc. Recommend surgical lubricants such as instrument milk are best. Do not use WD-40, oil or other industrial lubricants. Put instruments up for autoclaving either individually or in sets.

Inspection and Function Testing

Visually inspect to ensure all contamination has been removed. Check for distortion, damage and wear. Cutting edges should be free of defects. Discard damaged, worn or corroded instruments. Ensure that detachable tips are secure.

Drying

Dry using paper toweling or dry heat not exceeding 134 °C.

Packaging for Sterilization

Instruments may be loaded into dedicated instruments trays or general purpose sterilization trays, ensure that cutting edges are protected. For individual Instruments disposable paper or plastic pouches are ideal. Make sure you use a wide enough pouch (at least 10 cm) for instruments with ratchet locks such as needle holders and hemostats so the instruments can be sterilized in an open (unlocked) position.

Sterilization

The instruments can be sterilized in designated container systems or general purpose sterilization containers individually packed in the standard polyethylene or Sterilization must be done by steam with fractionated pre-vacuum. Based on DIN EN 285 and DIN EN 554 is

recommended at a temperature of 134 °C, a contact time of at least 5 minutes. Effective steam sterilization can be achieved with the following cycle: Temperature cycle type minimum exposure time Pre-vacuum to 134 °C - 137 °C for 5 minutes. To avoid damage the instrument, the temperature of 180 °C must not be exceeded.

Storage

Sterile, packaged instruments should be stored so that they are protected from dust, moisture, insects and vermin, and temperature extremes.

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