Instruction for use

This package insert is not a reference to surgical techniques. It is designed to assist in using this product.

Product description

Simplee Suture Material, absorbable is a synthetic absorbable sterile surgical suture composed of polyglactin, a copolymer made of 90% glycolide and 10% L-lactide. The coating (<1%) of the braided sutures is composed of poly-glycolid-co-L-lactide and calcium stearate. The braided sutures are dyed in violet. Sutures are also available in the undyed form. Progressive loss of tensile strength, mass and eventual absorption of the sutures occur by means of hydrolysis. The loss of effective tensile strength is approximately 40% - 50% 21 days post-implantation. The absorption is complete between 56 and 70 days.

Simplee Suture Material, absorbable is available in a range of gauge sizes and lengths, non-needled or attached to atraumatic stainless steel needles of varying types. The sutures are individually packed in sterile foil packs.

Simplee Suture Material, absorbable complies the requirements of the European Pharmacopoeia and the United States Pharmacopoeia.

Indications

Simplee Suture Material, absorbable is indicated for use in general soft tissue approximation and/or ligation.

The selection of the suture should be based on the patient's condition, the surgical experience, the surgical technique to be employed and the wound to be treated.

Contra-Indications

Simplee Suture Material, absorbable should not be used:

- in tissues which may undergo expansion, stretching or distension or which may require additional support,
- in cardiovascular and neurological tissues,
- in patients that are allergic to suture material.

The use of this suture can be improper in elders, patients that are not properly fed, weak patients, and patients under certain conditions that could delay the healing of the wound.

Adverse reactions

Adverse reactions associated with this device include: wound dehiscence, failure to provide adequate wound support in sites where expansion, stretching or distension occur or in patients with delayed wound healing, transient local irritation at the wound site, transient inflammatory foreign body response, induration, erythema at the wound side and calculi formation. Like all foreign bodies Simplee Suture Material, absorbable may potentiate an existing infection.

Handling instructions/ Warnings/ Precautions

Users should be familiar with surgical procedures and techniques and tie techniques involving absorbable sutures before employing Simplee Suture Material, absorbable suture for wound closure, as risk of wound dehiscence may vary with the site of application and the suture material used.

Particular factors of each patient should be considered regarding healing process in vivo.

Skin, vaginal and conjunctival sutures which must remain in place longer than 7 days may cause localised irritation and should be snipped off or removed as indicated. Place subcuticular sutures as deeply as possible to minimize the erythema and induration. Avoid prolonged contact of any suture with salt solutions (urinary or biliary), as it may result in calculus formation. Contaminated or infected wounds should be managed with acceptable surgical practice. A supplemental nonabsorbable suture or an immobilisation of joints by an external support should be considered by the surgeon in the closure of the sites which may require additional support (e.g. expansion, distension, etc.). In tissues with poor blood supply suture extrusion and delayed absorption may occur.

Open the outer bag of the foil pack by pulling its two loose ends away from each other – in such a way that the sterile bag containing the suture material falls out undamaged, onto a sterile surface. The sterile bag should then only be touched with sterile gloves and opened by tearing the upper part of the bag.

The product is intended for single use!

When suturing, remember:

- Always double check that you have tied a tight knot. Additional throws should be considered in some surgical circumstances.
- Users should exercise caution when handling surgical needles to avoid inadvertent needle stick injury, which may for example result in transmission of blood-borne pathogens.
- When suturing be careful to avoid damage (suture damage, needle damage, etc.) with surgical instruments (e.g. needle holder, forceps, etc.) or needle.
- Hold the needle in the area between the third and the half of the distance between the extreme of the union with the thread and the tip. If you hold it near to the tip, you could damage the functional integrity of the needle or break it. If you hold it by the union zone with the thread, the needle body could be broken. When the needle function is altered, it could generate a decrease in its folding resistance and break.
- Discard any unused suture material remaining. Suture needles should be discarded into separated "sharps" container. Dispose of contaminated devices and packaging materials utilizing standard hospital procedures and universal precautions for bio hazardous waste.
- Do not use rusty needles! Use of rusty needles may cause infection and other harms.
Do not resterilise!
Do not use if the package is opened or damaged.
This device was designed, tested and manufactured for single patient use only. Reuse, use of the device with opened or damaged packaging, reprocessing and/or resterilisation of this device may lead to its failure and subsequent injury, illness or death of the patient and/or create the risk of contamination and patient infection, illness or death of the patient.

Sterilisation

Simplee Suture Material, absorbable is sterilised by ethylene oxide.

Storage

The product should be protected from direct sunlight and heat, and stored in its original packaging in a clean, dry room at a temperature from 5°C to 30°C.

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LIFCO Dental AB · Verkstärestegatan 1