



SIMPLEE®

All you need. No more. But no less.

NEW

IN THE PRODUCT RANGE



SIMPLEE® NANO COMPOSITE

Light-curing nano composite, radiopaque, with ultra-fine glass fillers for extremely homogeneous and highly aesthetic restorations. The glass fillers enable optimal color matching of the composites.

Indication: Front and posterior tooth restorations of classes I-V (the color Universal is only suitable for the side tooth area as it is less translucent than the other colors). Inlays, onlays, veneers. Extended fissure sealing for molars and premolars. Core build-up. Splinting of loose teeth. Corrections of shape and color to improve aesthetics.

Farbe	Syringe, 4 g	Compules, 20 x 0,3 g
A2	156 130	156 136
A3	156 131	156 137
A3,5	156 132	156 138
B2	156 133	–
B3	156 134	–
Universal	156 135	156 139