



SIMPLEE®

All you need. No more. But no less.

NEW

IN THE PRODUCT RANGE



SIMPLEE® RECIPFILE V

Nickel-titanium files for mechanical root canal preparation with excellent cutting performance. Preparation by reciprocating movement. The files have very good flexibility and are therefore ideal for narrow, curved canals. The root canal can be optimally prepared with little effort using just a few instruments. Only one setting is required on the motor for the different file sizes: Rotation: 350 RPM; 3 – max. 4 Nm.

6 pieces

Size	R25	R40	R50
21 mm	156 640	156 643	156 646
25 mm	156 641	156 644	156 647
31 mm	156 642	156 645	156 648