



# SIMPLEE®

All you need. No more. But no less.

## NEW

## IN THE PRODUCT RANGE



### SIMPLEE® RECIPFILE VB

Blue nickel-titanium files with increased fracture resistance due to heat treatment of the root canal instruments. For mechanical root canal preparation with excellent cutting performance. Preparation by reciprocating movement. The files have very good flexibility and are therefore very suitable for narrow, curved canals. With just a few instruments, the root canal can be optimally prepared without great effort. Only one setting is required on the motor for the different file sizes: Rotation: 350 RPM; 3 – max.4 Nm.

#### 6 pieces

Size	R25	R40	R50
21 mm	<b>156 649</b>	<b>156 652</b>	<b>156 655</b>
25 mm	<b>156 650</b>	<b>156 653</b>	<b>156 656</b>
31 mm	<b>156 651</b>	<b>156 654</b>	<b>156 657</b>