



SIMPLEE®

All you need. No more. But no less.

NEW

IN THE PRODUCT RANGE



SIMPLEE® RETRACTION CORD

Unimpregnated, knitted retraction cord made of cotton. Due to the special knitting, small loops that form an interlocking chain, the threads can absorb a large amount of hemostatics and ensure good tissue displacement.

Thread length: 275 cm, container with integrated cutter in the lid.

Size	275 cm
------	--------

0	156 740
---	---------

00	156 741
----	---------

000	156 742
-----	---------

1	156 743
---	---------

2	156 744
---	---------