




SIMPLEE®

All you need. No more. But no less.

NEW

IN THE PRODUCT RANGE



SIMPLEE® W+ FILES

Nickel-titanium files for the mechanical, reciprocal preparation of the root canal. The W+ files are particularly flexible and also adapt perfectly to highly curved root canals. The files cover a wide range of canal morphologies. You can cover almost all root canal anatomies with just a few instruments. Only one setting is required on the motor for the different file sizes: Rotation: 350 RPM; 3 – max.4 Nm.

6 pieces

Size	large 45	medium 35	primary 25	small 20
21 mm	156 661	156 664	156 667	156 670
25 mm	156 662	156 665	156 668	156 671
31 mm	156 663	156 666	156 669	156 672